

Glypican 3 antibody



Catalog Number: 111089

Product name

Glypican 3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Glypican 3 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 66 kDa

Preparation

This antibody was obtained by immunization of Glypican 3 recombinant protein (Accession Number: NM_004484). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: N/A

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111089(GPC3 antibody) at dilution of 1:500