

# Glypican 4 antibody



Catalog Number: 111090

---

## Product name

Glypican 4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Glypican 4 Rabbit Polyclonal antibody. Positive WB detected in A549 cells, human brain tissue, mouse kidney tissue. Observed molecular weight by Western-blot: 62 kDa

## Preparation

This antibody was obtained by immunization of Glypican 4 recombinant protein (Accession Number: NM\_001448). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:111090(GPC4 antibody) at dilution of 1:200