

GMEB2 antibody



Catalog Number: 111097

Product name

GMEB2 antibody

Recommended Dilution:

WB: 1:200-1:2000

Immunogen

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GMEB2 Rabbit Polyclonal antibody. Positive IHC detected in human cervical cancer tissue. Positive WB detected in human placenta tissue, mouse ovary tissue, mouse thymus tissue. Observed molecular weight by Western-blot: 79 kDa or 67 kDa

Preparation

This antibody was obtained by immunization of GMEB2 recombinant protein (Accession Number: XM_047440106). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


ELISA, WB, IHC

Dilutions


Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:111097(GMEB2 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human cervical cancer using Catalog No:111097(GMEB2 antibody) at dilution of 1:100 (under 10x lens)



GMEB2 antibody



Catalog Number: 111097

Immunohistochemical of paraffin-embedded
human cervical cancer using Catalog

No:111097(GMEB2 antibody) at dilution of 1:100
(under 40x lens)