

GPR123-Specific antibody



Catalog Number: 111108

Product name

GPR123-Specific antibody

Specificity

Human; other species not tested.

Antibody description

GPR123-Specific Rabbit Polyclonal antibody.
Positive IHC detected in human testis tissue.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001083909).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:111108(GPR123-Specific Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:111108(GPR123-Specific Antibody) at dilution of 1:50 (under 40x lens)