

GPR132 antibody



Catalog Number: 111110

Product name

GPR132 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GPR132 Rabbit Polyclonal antibody. Positive WB detected in mouse spleen tissue. Positive IHC detected in human kidney tissue.

Preparation

This antibody was obtained by immunization of GPR132 recombinant protein (Accession Number: NM_001278694). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:111110(GPR132 Antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human kidney using Catalog No:111110(GPR132 antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:111110(GPR132 antibody) at dilution of 1:50 (under 10x lens)