

# GPR132 antibody

Catalog Number: 111110

## Product name

GPR132 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

GPR132 Rabbit Polyclonal antibody. Positive WB detected in mouse spleen tissue. Positive IHC detected in human kidney tissue.

## Preparation

This antibody was obtained by immunization of GPR132 recombinant protein (Accession Number: NM\_001278694). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


## Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:111110(GPR132 Antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human kidney using Catalog No:111110(GPR132 antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:111110(GPR132 antibody) at dilution of 1:50 (under 10x lens)