

GPR162 antibody



Catalog Number: 111117

Product name

Polyclonal

GPR162 antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

GPR162 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, HeLa cells, human brain tissue, Y79 cells. Observed molecular weight by Western-blot: 33 kDa, 64 kDa

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of GPR162 recombinant protein (Accession Number: NM_014449). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:111117(GPR162 antibody) at dilution of 1:500