

# GPR162 antibody



Catalog Number: 111117

---

## Product name

GPR162 antibody

## Specificity

Human; other species not tested.

## Antibody description

GPR162 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, HeLa cells, human brain tissue, Y79 cells. Observed molecular weight by Western-blot: 33 kDa, 64 kDa

## Preparation

This antibody was obtained by immunization of GPR162 recombinant protein (Accession Number: NM\_014449). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:111117(GPR162 antibody) at dilution of 1:500