

GPR21 antibody

Catalog Number: 111123

Product name

GPR21 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GPR21 Rabbit Polyclonal antibody. Positive IHC detected in human stomach tissue, human prostate hyperplasia tissue.

Preparation

This antibody was obtained by immunization of GPR21 recombinant protein (Accession Number: NM_005294). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:50-1:500

Validations



Immunohistochemistry of paraffin-embedded human stomach tissue slide using Catalog No:111123(GPR21 Antibody) at dilution of 1:100 (under 10x lens).



Immunohistochemistry of paraffin-embedded human stomach tissue slide using Catalog No:111123(GPR21 Antibody) at dilution of 1:100 (under 40x lens).