

# GPR37/Pael-R antibody



Catalog Number: 111127

## Product name

Polyclonal

GPR37/Pael-R antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

GPR37/Pael-R Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 50-55kd

## Dilutions

Recommended Dilution:

## Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of GPR37/Pael-R recombinant protein (Accession Number: NM\_005302). Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111127(GPR37 antibody) at dilution of 1:1600