

GPR58/TAAR2 antibody



Catalog Number: 111132

Product name

GPR58/TAAR2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GPR58/TAAR2 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue, mouse brain tissue, mouse cerebellum tissue. Observed molecular weight by Western-blot: 42 kDa

Preparation

This antibody was obtained by immunization of GPR58/TAAR2 recombinant protein (Accession Number: NM_014626). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111132(TAAR2 Antibody) at dilution of 1:600