

GPR58/TAAR2 antibody



Catalog Number: 111132

Product name

Polyclonal

GPR58/TAAR2 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

GPR58/TAAR2 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue, mouse brain tissue, mouse cerebellum tissue. Observed molecular weight by Western-blot: 42 kDa

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of GPR58/TAAR2 recombinant protein (Accession Number: NM_014626). Purification method: Antigen affinity purified.

Formulation


PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111132(TAAR2 Antibody) at dilution of 1:600