

Product name

GPRC5D antibody

Specificity

Human; other species not tested.

Antibody description

GPRC5D Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue.

Preparation

This antibody was obtained by immunization of GPRC5D recombinant protein (Accession Number: NM_018654). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC


Dilutions

Recommended Dilution:

IHC: 1:50-1:500

Validations

Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:111140(GPRC5D Antibody) at dilution of 1:100 (under 10x lens).



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:111140(GPRC5D Antibody) at dilution of 1:100 (under 40x lens).