

# GRAP2 antibody



Catalog Number: 111162

## Product name

GRAP2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

GRAP2 Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue. Positive WB detected in Transfected HEK-293 cells.

## Preparation

This antibody was obtained by immunization of GRAP2 recombinant protein (Accession Number: NM\_001291824). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:111162 (GRAP2 antibody) at dilution of 1:50 (under 10x lens)



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111162 (GRAP2 Antibody) at dilution of 1:1000