

GSTO2 antibody



Catalog Number: 111190

Product name

GSTO2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GSTO2 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, HeLa cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 28 kDa

Preparation

This antibody was obtained by immunization of GSTO2 recombinant protein (Accession Number: XM_011539270). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:111190(GSTO2 antibody) at dilution of 1:300



Immunofluorescent analysis of (-20oc Ethanol) fixed HeLa cells using Catalog No:111190(GSTO2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)