

# GSTT2 antibody



Catalog Number: 111193

## Product name

GSTT2 antibody

## Immunogen

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

GSTT2 Rabbit Polyclonal antibody. Positive IF detected in L02 cells. Positive WB detected in rat liver tissue, L02 cells, mouse liver tissue, mouse lung tissue. Observed molecular weight by Western-blot: 28 kDa

## Preparation

This antibody was obtained by immunization of GSTT2 recombinant protein (Accession Number: NM\_001080843). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

## Validations



rat liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111193(GSTT2 Antibody) at dilution of 1:1500



Immunofluorescent analysis of L02 cells using Catalog No:111193(GSTT2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)