

GRK3 antibody

Catalog Number: 111213

Product name

GRK3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GRK3 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells. Observed molecular weight by Western-blot: 70kd and 90kd

Preparation

This antibody was obtained by immunization of GRK3 recombinant protein (Accession Number: XM_047441167). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations

K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:111213(ADRBK2 antibody) at dilution of 1:400