

GRK5 antibody



Catalog Number: 111214

Product name

Polyclonal

GRK5 antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

GRK5 Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue. Observed molecular weight by Western-blot: 68kd

Dilutions

Recommended Dilution:

Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of GRK5 recombinant protein (Accession Number: NM_005308). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:111214(GRK5 antibody) at dilution of 1:400