

GTF2E1 antibody



Catalog Number: 111229

Product name

GTF2E1 antibody

Specificity

Human; other species not tested.

Antibody description

GTF2E1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 50 kDa

Preparation

This antibody was obtained by immunization of GTF2E1 recombinant protein (Accession Number: NM_005513). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111229(GTF2E1 Antibody) at dilution of 1:3000