

GTPBP1 antibody

Catalog Number: 111241

Product name

GTPBP1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GTPBP1 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, Jurkat cells. Positive IP detected in mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 72 kDa

Preparation

This antibody was obtained by immunization of GTPBP1 recombinant protein (Accession Number: NM_004286). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:111241(GTPBP1 antibody) at dilution of 1:600



IP Result of anti-GTPBP1 (IP:Catalog No:111241, 4ug; Detection:Catalog No:111241 1:500) with mouse skeletal muscle tissue lysate 4000ug.