

GUCY1B3 antibody



Catalog Number: 111249

Product name

GUCY1B3 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

GUCY1B3 Rabbit Polyclonal antibody. Positive IF detected in Hela. Positive WB detected in human placenta tissue. Observed molecular weight by Western-blot: 71kd

Preparation

This antibody was obtained by immunization of GUCY1B3 recombinant protein (Accession Number: NM_000857). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:111249(GUCY1B3 antibody) at dilution of 1:500



Immunofluorescent analysis of Hela cells, using GUCY1B3 antibody Catalog No:111249 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).