

# Harmonin antibody



Catalog Number: 111254

## Product name

Harmonin antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

Harmonin Rabbit Polyclonal antibody. Positive WB detected in human small intestine tissue. Positive IHC detected in human pancreas cancer tissue. Observed molecular weight by Western-blot: 73kd

## Preparation

This antibody was obtained by immunization of Harmonin recombinant protein (Accession Number: NM\_001297764). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



human small intestine tissue were subjected to SDS PAGE followed by western blot with Catalog No:111254(USH1C antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:111254(USH1C antibody) at dilution of 1:50 (under 10x lens)