

# HAS3 antibody



Catalog Number: 111259

## Product name

HAS3 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

HAS3 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, mouse liver tissue. Observed molecular weight by Western-blot: 63kd

## Preparation

This antibody was obtained by immunization of HAS3 recombinant protein (Accession Number: NM\_138612). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:111259(HAS3 antibody) at dilution of 1:1000