

# HBA1-Specific antibody



Catalog Number: 111265

## Product name

Polyclonal

HBA1-Specific antibody

## Ig Type

## Specificity

Rabbit IgG

Human; other species not tested.

## Applications

## Antibody description

ELISA, WB

HBA1-Specific Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, human spleen tissue. Observed molecular weight by Western-blot: 14-16kd

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_000558). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:111265(HBA1-Specific antibody) at dilution of 1:300