

# HBQ1-Specific antibody



Catalog Number: 111273

## Product name

HBQ1-Specific antibody

## Specificity

Human; other species not tested.

## Antibody description

HBQ1-Specific Rabbit Polyclonal antibody.  
Positive WB detected in recombinant protein.  
Observed molecular weight by Western-blot: 42kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_005331).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



WB results of Catalog No:111273 (HBQ1 antibody) at dilution of 1:2000.