

# HBXIP antibody

Catalog Number: 111275

## Product name

HBXIP antibody

WB: 1:200-1:2000

IHC: 1:20-1:200

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

HBXIP Rabbit Polyclonal antibody. Positive FC detected in MCF-7 cells. Positive IHC detected in human hepat°Cirrhosis tissue, human breast cancer tissue, human heart tissue, human liver cancer tissue. Positive WB detected in Transfected HEK-293 cells. Observed molecular weight by Western-blot: 22 kDa

## Preparation

This antibody was obtained by immunization of HBXIP recombinant protein (Accession Number: NM\_006402). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IHC, FC


## Dilutions

Recommended Dilution:


## Validations



Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog No:111275(HBXIP antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog No:111275(HBXIP antibody) at dilution of 1:50 (under 40x lens)



# HBXIP antibody



Catalog Number: 111275

---

Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111275(HBXIP Antibody) at dilution of 1:800

1X10<sup>6</sup> MCF-7 cells were stained with 0.2ug HBXIP antibody (Catalog No:111275, red) and control antibody (blue). Fixed with 4% PFA blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.

