

# HBZ-Specific antibody

Catalog Number: 111277

## Product name

HBZ-Specific antibody

## Specificity

Human; other species not tested.

## Antibody description

HBZ-Specific Rabbit Polyclonal antibody. Positive WB detected in recombinant protein. Observed molecular weight by Western-blot: 42kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_001003938). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



WB results of Catalog No:111277 (HBZ antibody) at dilution of 1:2000.