

# Hephaestin antibody



Catalog Number: 111292

---

## Product name

Hephaestin antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Hephaestin Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue.

## Preparation

This antibody was obtained by immunization of Hephaestin recombinant protein (Accession Number: NM\_001367233). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions

Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:111292(HEPH antibody) at dilution of 1:50 (under 10x lens)