

HERC5 antibody

Catalog Number: 111297

Product name

HERC5 antibody

Specificity

Human; other species not tested.

Antibody description

HERC5 Rabbit Polyclonal antibody. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 117 kDa

Preparation

This antibody was obtained by immunization of HERC5 recombinant protein (Accession Number: NM_016323). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IP

Dilutions

Recommended Dilution:

IP: 1:200-1:1000

Validations



IP Result of anti-HERC5 (IP:Catalog No:111297, 4ug; Detection:Catalog No:111297 1:300) with HEK-293 cells lysate 2800ug.