

# HES5 antibody

Catalog Number: 111300

## Product name

HES5 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

HES5 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue. Positive IHC detected in mouse brain tissue. Observed molecular weight by Western-blot: 41 kDa

## Preparation

This antibody was obtained by immunization of HES5 recombinant protein (Accession Number: BC087840). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

## Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111300(HES5 Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:111300(HES5 Antibody) at dilution of 1:50 (under 40x lens)