

# HESX1 antibody

Catalog Number: 111301

---

**Product name**

HESX1 antibody

**Specificity**

Human, Mouse, Rat; other species not tested.

**Antibody description**

HESX1 Rabbit Polyclonal antibody. Positive IHC detected in human brain.

**Preparation**

This antibody was obtained by immunization of HESX1 recombinant protein (Accession Number: XM\_005265526). Purification method: Antigen affinity purified.

**Formulation**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

**Storage**

Store at -20°C. DO NOT ALIQUOT

**Clonality**

Polyclonal

**Ig Type**

Rabbit IgG


**Applications**

ELISA, IHC

**Dilutions**

Recommended Dilution:

IHC: 1:20-1:200

**Validations**

Immunohistochemical of paraffin-embedded human brain using Catalog No:111301(HESX1 antibody) at dilution of 1:50 (under 10x lens)