

# GYG2 antibody



Catalog Number: 111319

---

## Product name

GYG2 antibody

## Specificity

human; other species not tested.

## Antibody description

GYG2 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, HEK-293 cells, human liver tissue. Observed molecular weight by Western-blot: 45 kDa, 66-70kDa

## Preparation

This antibody was obtained by immunization of GYG2 recombinant protein (Accession Number: XM\_011545600). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:111319(GYG2 antibody) at dilution of 1:500