

HIC2 antibody

Catalog Number: 111343

Product name

HIC2 antibody

Specificity

Human; other species not tested.

Antibody description

HIC2 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Positive IP detected in Jurkat cells. Observed molecular weight by Western-blot: 66 kDa

Preparation

This antibody was obtained by immunization of HIC2 recombinant protein (Accession Number: XM_011530008). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111343(HIC2 Antibody) at dilution of 1:300



IP Result of anti-HIC2 (IP:Catalog No:111343, 4ug; Detection:Catalog No:111343 1:300) with Jurkat cells lysate 2400ug.