

# HIST1H1T antibody



Catalog Number: 111359

## Product name

HIST1H1T antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

HIST1H1T Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, mouse kidney tissue, mouse testis tissue. Observed molecular weight by Western-blot: 30 kDa

## Preparation

This antibody was obtained by immunization of HIST1H1T recombinant protein (Accession Number: NM\_005323). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions


Recommended Dilution:

WB: 1:200-1:1000

## Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:111359(HIST1H1T Antibody) at dilution of



1:300

mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:111359(HIST1H1T Antibody) at dilution of 1:300