

HIST1H4F antibody



Catalog Number: 111362

Product name

HIST1H4F antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HIST1H4F Rabbit Polyclonal antibody. Positive WB detected in human spleen tissue. Observed molecular weight by Western-blot: 11 kDa

Preparation

This antibody was obtained by immunization of HIST1H4F recombinant protein (Accession Number: NM_003540). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



human spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:111362(HIST1H4F antibody) at dilution of 1:300