

# HIST2H2AA4 antibody



Catalog Number: 111363

## Product name

HIST2H2AA4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

HIST2H2AA4 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, Jurkat cells, K-562 cells, mouse liver tissue, rat liver tissue. Positive IHC detected in human liver cancer tissue. Observed molecular weight by Western-blot: 14 kDa

## Preparation

This antibody was obtained by immunization of HIST2H2AA4 recombinant protein (Accession Number: BC001629). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:111363(HIST2H2AA4 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:111363(HIST2H2AA4 Antibody) at dilution of 1:50 (under 40x lens)