

HDAC4 antibody

Catalog Number: 111377

Product name

HDAC4 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HDAC4 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, Jurkat cells. Positive IP detected in HeLa cells. Positive IHC detected in human kidney tissue, human cervical cancer tissue. Observed molecular weight by Western-blot: 119-140kd

Preparation

This antibody was obtained by immunization of HDAC4 recombinant protein (Accession Number: BC039904). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111377(HDAC4 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human kidney using Catalog No:111377(HDAC4 antibody) at dilution of 1:100 (under 40x lens)




IP Result of anti-HDAC4 (IP:Catalog No:111377, 4ug; Detection:Catalog No:111377 1:500) with

HDAC4 antibody

Catalog Number: 111377

HeLa cells lysate 4000ug.



Immunohistochemistry of paraffin-embedded human cervical cancer tissue slide using Catalog No:111377(HDAC4 Antibody) at dilution of 1:50 (under 40x lens)