

# HDAC4 antibody



Catalog Number: 111377

## Product name

HDAC4 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

HDAC4 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, Jurkat cells. Positive IP detected in HeLa cells. Positive IHC detected in human kidney tissue, human cervical cancer tissue. Observed molecular weight by Western-blot: 119-140kd

## Preparation

This antibody was obtained by immunization of HDAC4 recombinant protein (Accession Number: BC039904). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000


## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111377(HDAC4 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human kidney using Catalog No:111377(HDAC4 antibody) at dilution of 1:100 (under 40x lens)



IP Result of anti-HDAC4 (IP:Catalog No:111377, 4ug; Detection:Catalog No:111377 1:500) with


# HDAC4 antibody



Catalog Number: 111377

---

HeLa cells lysate 4000ug.



Immunohistochemistry of paraffin-embedded human cervical cancer tissue slide using Catalog No:111377(HDAC4 Antibody) at dilution of 1:50 (under 40x lens)