

# HEATR4 antibody



Catalog Number: 111393

## Product name

HEATR4 antibody

## Immunogen

[Human HEATR4 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

HEATR4 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue, human breast cancer tissue.

## Preparation

This antibody was obtained by immunization of HEATR4 recombinant protein (Accession Number: BC047590). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human testis slide using Catalog No:111393(HEATR4 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human testis slide using Catalog No:111393(HEATR4 Antibody) at dilution of 1:50 (under 40x lens)