

HELT antibody



Catalog Number: 111401

Product name

HELT antibody

Immunogen

[Human HELT Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HELT Rabbit Polyclonal antibody. Positive WB detected in rat heart tissue, mouse heart tissue. Observed molecular weight by Western-blot: 27 kDa

Preparation

This antibody was obtained by immunization of HELT recombinant protein (Accession Number: NM_001300782). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



rat heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:111401(HELT Antibody) at dilution of 1:1000