

HIVEP1 antibody

Catalog Number: 111410

Product name

HIVEP1 antibody

Specificity

Human; other species not tested.

Antibody description

HIVEP1 Rabbit Polyclonal antibody. Positive IHC detected in human stomach tissue, human testis tissue.

Preparation

This antibody was obtained by immunization of HIVEP1 recombinant protein (Accession Number: XM_047418700). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human stomach tissue slide using Catalog No:111410(HIVEP1 Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:111410(HIVEP1 Antibody) at dilution of 1:50 (under 40x lens)