

HLA-DQB2 antibody



Catalog Number: 111424

Product name

HLA-DQB2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

HLA-DQB2 Rabbit Polyclonal antibody. Positive WB detected in Raji cells, HepG2 cells, K-562 cells. Observed molecular weight by Western-blot: 27kd,31kd

Preparation

This antibody was obtained by immunization of HLA-DQB2 recombinant protein (Accession Number: BC031995). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:111424(HLA-DQB2 antibody) at dilution of 1:500