

# HLA-DRA antibody



Catalog Number: 111425

## Product name

HLA-DRA antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

HLA-DRA Rabbit Polyclonal antibody. Positive IHC detected in human tonsil.

## Preparation

This antibody was obtained by immunization of HLA-DRA recombinant protein (Accession Number: NM\_019111). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human tonsil using Catalog No:111425(HLA-DRA antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human tonsil using Catalog No:111425(HLA-DRA antibody) at dilution of 1:50 (under 40x lens)