

# HLA-DRB5 antibody



Catalog Number: 111428

## Product name

HLA-DRB5 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

HLA-DRB5 Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue. Positive WB detected in HepG2 cells, HeLa cells, HL-60 cells. Observed molecular weight by Western-blot: 30 kDa

## Preparation

This antibody was obtained by immunization of HLA-DRB5 recombinant protein (Accession Number: NM\_002125). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:111428(HLA-DRB5 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human tonsillitis using Catalog No:111428(HLA-DRB5 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human tonsillitis using Catalog No:111428(HLA-DRB5 antibody) at dilution of 1:50 (under 40x

# HLA-DRB5 antibody



Catalog Number: 111428

---

lens)