

HLA-E antibody



Catalog Number: 111429

Product name

HLA-E antibody

Specificity

Human; other species not tested.

Antibody description

HLA-E Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, human kidney tissue. Observed molecular weight by Western-blot: 40kd

Preparation

This antibody was obtained by immunization of HLA-E recombinant protein (Accession Number: NM_005516). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111429(HLA-E antibody) at dilution of 1:500