

HLCS antibody



Catalog Number: 111432

Product name

HLCS antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HLCS Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, human kidney tissue, human liver tissue. Positive IF detected in HEK-293 cells. Observed molecular weight by Western-blot: 76 kDa

Preparation

This antibody was obtained by immunization of HLCS recombinant protein (Accession Number: XM_047440754). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111432(HLCS antibody) at dilution of 1:700



Immunofluorescent analysis of HEK-293 cells using Catalog No:111432(HLCS Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)