

HLF antibody



Catalog Number: 111433

Product name

HLF antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

HLF Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, HEK-293 cells, HepG2 cells, K-562 cells, mouse brain tissue. Observed molecular weight by Western-blot: 66kd

Preparation

This antibody was obtained by immunization of HLF recombinant protein (Accession Number: NM_002126). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111433(HLF antibody) at dilution of 1:800