

HLTF antibody

Catalog Number: 111434

Product name

HLTF antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HLTF Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, K-562 cells. Observed molecular weight by Western-blot: 114 kDa

Preparation

This antibody was obtained by immunization of HLTF recombinant protein (Accession Number: XM_017007079). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111434(HLTF antibody) at dilution of 1:300



IP Result of anti-HLTF (IP:Catalog No:111434, 4ug; Detection:Catalog No:111434 1:500) with HeLa cells lysate 2400ug.