

# HMBOX1 antibody

Catalog Number: 111436

## Product name

HMBOX1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

HMBOX1 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 47 kDa

## Preparation

This antibody was obtained by immunization of HMBOX1 recombinant protein (Accession Number: BC009259). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111436(HMBOX1 antibody) at dilution of



1:1000

IP Result of anti-HMBOX1 (IP:Catalog No:111436, 4ug; Detection:Catalog No:111436 1:1000) with mouse brain tissue lysate 3440ug.