

# HOXD12 antibody



Catalog Number: 111437

## Product name

Polyclonal

HOXD12 antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse; other species not tested.

## Applications

## Antibody description

ELISA, WB

HOXD12 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 30 kDa

## Dilutions

Recommended Dilution:


## Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of HOXD12 recombinant protein (Accession Number: BC121103). Purification method: Antigen affinity purified.

## Validations

## Formulation



PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111437(HOXD12 Antibody) at dilution of 1:600

## Clonality

Store at -20°C. DO NOT ALIQUOT